

CIRS Tissue Equivalent Materials

Model XXX



FOR DOSIMETRY & IMAGE QUALITY STUDIES

CIRS Tissue Simulation Technology allows users to study the interaction between ionizing radiation and tissue with safety and precision. Our tissue equivalent materials (TEMs) have been validated through independent studies, continuous monitoring and worldwide use for over 35 years. All TEMs produced by CIRS meet International Commission on Radiation Units and Measurements (ICRU) Report 44 standards, which require that TEMs match the linear attenuation and stopping power coefficient of their reference tissues to within +/- 1%. (Physical densities may vary by more than 1%, as it is sometimes necessary to adjust the physical density to ensure that the linear attenuation coefficient and electron density meets ICRU Report 44 requirements.)

CIRS Tissue Equivalent Materials are available in sizes ranging from 10 x 10 cm to 30 cm x 30 cm and thicknesses of 0.1 cm through 5 cm. Materials in different shapes and sizes are also available upon request.

CIRS TEM's are easily machined and can be glued together to create a thicker bolus of material. Other formulations are available upon request.

Benefits

- Meets ICRU Report 44 standards for tissue equivalence in radiation therapy and diagnostic imaging
- Applies to photon, electron and proton beams
- Covers therapeutic and diagnostic energy ranges (see page 2)
- Tissue types available:
 - Lung (Inhale, exhale and medium)
 - Bone (Cortical, trabecular and average)
 - Breast (Adipose, glandular and mixed)
 - Organ Soft Tissue (Adipose Av. soft tissue, bladder, blood, brain white Mt, brain gray Mt, brain average, muscle, kidney, pancreas, heart, liver, prostate, spleen, intestine)

深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机)

地址:深圳市福田区八卦二路八卦岭工业区615栋419

邮箱:hongqi@thingstet.com

网址:www.thingstest.com

CIRS

Tissue Simulation & Phantom Technology

TISSUE EQUIVALENT MATERIALS

DOSIMETRY


Slabs can be manufactured to accommodate a variety of detectors including ion chambers, diodes, MOSFETs, TLD and OSL dosimeters in standard or custom locations.

Standard detector location is in the middle section of the 2 cm slab in the geometrical center of the slab's area.

Location is scribed with a 10x10 cm field and external laser marks for easy positioning of slab against radiation beam's light field and lasers.



Photon Energy Coverage - TEMs



SPECIFICATIONS

TISSUE EQUIVALENT MATERIALS INCLUDE

| QTY | DESCRIPTION |
|-----|-------------------------------|
| 1 | Tissue Equivalent Slab or Set |
| 1 | User Guide |
| - | 60 Month Warranty |

See pages 3-4 for detailed technical specifications on each tissue type available.

深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机)

地址:深圳市福田区八卦二路八卦岭工业区615栋419

邮箱:hongqi@thingstet.com

网址:www.thingstest.com

TISSUE EQUIVALENT MATERIALS

LUNG

Lung - available in Inhale (0.2 g/cc), Medium (0.3 g/cc) or Exhale (0.5 g/cc).

| SLAB SIZE (CM) | THICKNESS (CM) | | | | |
|----------------|----------------|-----|-----|-----|-----|
| | 0.5 | 1.0 | 2.0 | 3.0 | 4.0 |
| 20 X 20 | | X | X | X | X |
| 30 X 30 | | X | X | X | X |

Technical data for the CIRS Lung-equivalent materials.

Recalculated Linear attenuation coefficients (cm^{-1}).

| En, MeV | Lung tissue (inhale) | | | Lung tissue (medium) | | | Lung tissue (exhale) | | |
|--|------------------------|--------|----------|------------------------|--------|----------|------------------------|--------|----------|
| | Reference ² | CIRS | Ratio, % | Reference ² | CIRS | Ratio, % | Reference ² | CIRS | Ratio, % |
| 0.04 | 0.0537 | 0.0524 | 97.6 | 0.0805 | 0.0785 | 97.5 | 0.1342 | 0.1313 | 97.8 |
| 0.06 | 0.0410 | 0.0411 | 100.2 | 0.0615 | 0.0613 | 99.7 | 0.1025 | 0.1012 | 98.7 |
| 0.08 | 0.0365 | 0.0367 | 100.6 | 0.0547 | 0.0550 | 100.5 | 0.0912 | 0.0904 | 99.1 |
| 0.10 | 0.0339 | 0.0341 | 100.6 | 0.0509 | 0.0513 | 100.8 | 0.0848 | 0.0840 | 99.1 |
| 0.20 | 0.0272 | 0.0274 | 100.7 | 0.0408 | 0.0413 | 101.2 | 0.0680 | 0.0675 | 99.3 |
| 0.40 | 0.0211 | 0.0212 | 100.5 | 0.0316 | 0.0320 | 101.3 | 0.0526 | 0.0523 | 99.4 |
| 0.60 | 0.0178 | 0.0179 | 100.6 | 0.0267 | 0.0270 | 101.1 | 0.0444 | 0.0442 | 99.5 |
| 0.80 | 0.0156 | 0.0157 | 100.6 | 0.0234 | 0.0237 | 101.3 | 0.0390 | 0.0388 | 99.5 |
| 1.00 | 0.0140 | 0.0141 | 100.7 | 0.0210 | 0.0213 | 101.4 | 0.0351 | 0.0349 | 99.4 |
| 2.00 | 0.0098 | 0.0099 | 101.0 | 0.0147 | 0.0149 | 101.4 | 0.0245 | 0.0243 | 99.2 |
| 4.00 | 0.0068 | 0.0068 | 100.0 | 0.0101 | 0.0102 | 101.0 | 0.0169 | 0.0167 | 98.8 |
| 6.00 | 0.0055 | 0.0055 | 100.0 | 0.0082 | 0.0083 | 101.2 | 0.0137 | 0.0135 | 98.5 |
| 8.00 | 0.0048 | 0.0048 | 100.0 | 0.0072 | 0.0072 | 100.0 | 0.0120 | 0.0117 | 97.5 |
| 10.0 | 0.0044 | 0.0043 | 97.7 | 0.0066 | 0.0066 | 100.0 | 0.0110 | 0.0106 | 96.4 |
| 20.0 | 0.0036 | 0.0035 | 97.2 | 0.0054 | 0.0053 | 98.1 | 0.0090 | 0.0085 | 94.4 |
| 30.0 | 0.0034 | 0.0032 | 94.1 | 0.0051 | 0.0050 | 98.0 | 0.0085 | 0.0079 | 92.9 |
| El. density, $\times 10^{23}, \text{cm}^{-3}$ | 0.663 | 0.668 | 100.8 | 0.995 | 1.008 | 101.3 | 1.658 | 1.648 | 99.4 |
| Density, g/cm^3 | 0.20 | 0.205 | | 0.30 | 0.31 | | 0.50 | 0.50 | |
| ORDER CODE | LAA | | | LG3 | | | LH | | |

1. ICRP 23, Report of the Task Group on Reference Man (1975).

2. Woodard, H.Q., White, D.R., The Composition of Body Tissues, The British Journal of Radiology (1986) 59: 1209-1219.

BONE

Bone - available in Average, Cortical, or Trabecular

| SLAB SIZE (CM) | THICKNESS (CM) | | | |
|----------------|----------------|-----|-----|-----|
| | 0.5 | 1.0 | 2.0 | 3.0 |
| 10 X 10 | | | | |
| 20 X 20 | X | X | X | X |
| 30 X 30 | X | X | X | X |

Technical data for the CIRS Bone-equivalent materials.

Recalculated linear attenuation coefficients (cm^{-1}).

| En, MeV | Spongiosa (trabecular) bone tissue | | | Average bone tissue* | | | Cortical bone tissue | | |
|--|------------------------------------|--------|----------|------------------------|--------|----------|------------------------|--------|----------|
| | Reference ² | CIRS | Ratio, % | Reference ² | CIRS | Ratio, % | Reference ² | CIRS | Ratio, % |
| 0.04 | 0.4546 | 0.4536 | 99.8 | 0.7884 | 0.7891 | 100.09 | 1.2783 | 1.2693 | 99.3 |
| 0.06 | 0.2802 | 0.2806 | 100.1 | 0.4244 | 0.4245 | 100.02 | 0.6046 | 0.6025 | 99.7 |
| 0.08 | 0.2296 | 0.2303 | 100.3 | 0.3251 | 0.3250 | 99.97 | 0.4282 | 0.4273 | 99.8 |
| 0.10 | 0.2058 | 0.2065 | 100.3 | 0.2822 | 0.2821 | 99.96 | 0.3561 | 0.3560 | 100.0 |
| 0.20 | 0.1588 | 0.1596 | 100.5 | 0.2098 | 0.2096 | 99.90 | 0.2517 | 0.2513 | 99.8 |
| 0.40 | 0.1223 | 0.1229 | 100.5 | 0.1605 | 0.1603 | 99.88 | 0.1903 | 0.1903 | 100.0 |
| 0.60 | 0.1031 | 0.1036 | 100.5 | 0.1351 | 0.1350 | 99.93 | 0.1600 | 0.1601 | 100.1 |
| 0.80 | 0.0905 | 0.0909 | 100.4 | 0.1186 | 0.1185 | 99.92 | 0.1403 | 0.1404 | 100.1 |
| 1.00 | 0.0813 | 0.0817 | 100.5 | 0.1066 | 0.1065 | 99.91 | 0.1260 | 0.1261 | 100.1 |
| 2.00 | 0.0568 | 0.0571 | 100.5 | 0.0746 | 0.0746 | 100.00 | 0.0884 | 0.0885 | 100.1 |
| 4.00 | 0.0393 | 0.0395 | 100.5 | 0.0521 | 0.0520 | 99.81 | 0.0626 | 0.0624 | 99.7 |
| 6.00 | 0.0322 | 0.0323 | 100.3 | 0.0431 | 0.0431 | 100.00 | 0.0525 | 0.0523 | 99.6 |
| 8.00 | 0.0284 | 0.0284 | 100.0 | 0.0383 | 0.0383 | 100.00 | 0.0473 | 0.0471 | 99.6 |
| 10.0 | 0.0260 | 0.0260 | 100.0 | 0.0355 | 0.0355 | 100.00 | 0.0444 | 0.0441 | 99.3 |
| 20.0 | 0.0216 | 0.0215 | 99.5 | 0.0305 | 0.0305 | 100.00 | 0.0397 | 0.0391 | 98.5 |
| 30.0 | 0.0206 | 0.0205 | 99.5 | 0.0296 | 0.0296 | 100.00 | 0.0394 | 0.0387 | 98.2 |
| El. density, $\times 10^{23}, \text{cm}^{-3}$ | 3.844 | 3.863 | 100.5 | 5.035 | 5.030 | 99.90 | 5.952 | 5.956 | 100.1 |
| Density, g/cm^3 | 1.18 | 1.20 | | 1.577 | 1.60 | | 1.92 | 1.93 | |
| ORDER CODE | SB-DTB | | | AB-IBN | | | CB-19F2 | | |

* The elemental composition and density were calculated from the average skeleton data minus cartilage.

1. ICRP 23, Report of the Task Group on Reference Man (1975).

2. Woodard, H.Q., White, D.R., The Composition of Body Tissues, The British Journal of Radiology (1986) 59: 1209-1219.

深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机)

地址:深圳市福田区八卦二路八卦岭工业区615栋419

邮箱:hongqi@thingstest.com

网址:www.thingstest.com

TISSUE EQUIVALENT MATERIALS

BREAST SOFT TISSUES

Adipose and Glandular Compositions - available in:

| | |
|----------------|------------------------------|
| 100% Adipose | 30% Glandular/ 70% Adipose |
| 100% Glandular | 50% Glandular / 50% Adipose |
| BR12 | 70% Glandular/ 30% Adipose |
| Swirled | 100% Glandular/ 100% Adipose |

| SLAB SIZE (CM) | THICKNESS (CM) | | | | |
|----------------|----------------|-----|-----|-----|-----|
| | 0.5 | 1.0 | 2.0 | 3.0 | 4.0 |
| 10 X 10 | X | X | X | X | X |
| 10 X 12.5 | X | X | X | | |
| 18 X 24 | | X | X | X | X |
| 24 X 30 | X | X | X | X | X |
| 30 X 30 | X | X | X | X | X |

Technical data for the CIRS Inc. breast tissue-equivalent materials.

Linear attenuation coefficients (cm^{-1}).

February 11, 2010

| En, MeV | Breast tissue - Glandular | | | Breast tissue- 70/30 (Glandular /Adipose) | | | Breast tissue - 50/50 (Glandular /Adipose) | | |
|--|---------------------------|--------|----------|---|--------|----------|--|--------|----------|
| | Reference ¹ | CIRS | Ratio, % | Reference ¹ | CIRS | Ratio, % | Reference ¹ | CIRS | Ratio, % |
| 0.01 | 4.9195 | 4.6315 | 94.15 | 4.2397 | 4.0281 | 95.01 | 3.8120 | 3.6543 | 95.86 |
| 0.02 | 0.7680 | 0.7680 | 100.00 | 0.6815 | 0.6802 | 99.81 | 0.6272 | 0.6262 | 99.84 |
| 0.04 | 0.2688 | 0.2708 | 100.74 | 0.2528 | 0.2536 | 100.32 | 0.2428 | 0.2433 | 100.21 |
| 0.06 | 0.2099 | 0.2103 | 100.19 | 0.2010 | 0.2008 | 99.90 | 0.1954 | 0.1952 | 99.90 |
| 0.08 | 0.1883 | 0.1883 | 100.00 | 0.1813 | 0.1808 | 99.72 | 0.1770 | 0.1765 | 99.72 |
| 0.10 | 0.1754 | 0.1753 | 99.94 | 0.1693 | 0.1687 | 99.65 | 0.1655 | 0.1650 | 99.70 |
| 0.20 | 0.1411 | 0.1410 | 99.93 | 0.1366 | 0.1360 | 99.56 | 0.1337 | 0.1333 | 99.70 |
| 0.40 | 0.1094 | 0.1093 | 99.91 | 0.1059 | 0.1055 | 99.62 | 0.1037 | 0.1034 | 99.71 |
| 0.60 | 0.0923 | 0.0923 | 100.00 | 0.0894 | 0.0891 | 99.66 | 0.0875 | 0.0873 | 99.77 |
| 0.80 | 0.0811 | 0.0810 | 99.88 | 0.0785 | 0.0782 | 99.62 | 0.0769 | 0.0767 | 99.74 |
| 1.00 | 0.0729 | 0.0728 | 99.86 | 0.0706 | 0.0703 | 99.58 | 0.0691 | 0.0689 | 99.71 |
| 2.00 | 0.0509 | 0.0508 | 99.80 | 0.0493 | 0.0491 | 99.59 | 0.0482 | 0.0481 | 99.79 |
| 4.00 | 0.0350 | 0.0348 | 99.43 | 0.0338 | 0.0336 | 99.41 | 0.0331 | 0.0329 | 99.40 |
| 6.00 | 0.0285 | 0.0281 | 98.60 | 0.0274 | 0.0271 | 98.91 | 0.0268 | 0.0266 | 99.25 |
| 8.00 | 0.0249 | 0.0245 | 98.39 | 0.0239 | 0.0236 | 98.74 | 0.0234 | 0.0231 | 98.72 |
| 10.0 | 0.0227 | 0.0222 | 97.80 | 0.0218 | 0.0214 | 98.17 | 0.0212 | 0.0209 | 98.58 |
| 20.0 | 0.0184 | 0.0177 | 96.20 | 0.0176 | 0.0170 | 96.59 | 0.0170 | 0.0166 | 97.65 |
| El. density, * 10^{23} , cm^{-3} | 3.445 | 3.443 | 99.94 | 3.336 | 3.325 | 99.67 | 3.267 | 3.258 | 99.72 |
| Density, g cm^{-3} | 1.04 | 1.048 | | 1.004 | 1.011 | | 0.982 | 0.99 | |

ORGAN SOFT TISSUES ICRU 46

Adipose Av. Soft Tissue, Bladder, Blood
 Brain White Mt, Brain Gray Mt, Brain Average
 Muscle, Kidney, Pancreas, Heart
 Liver, Prostate, Spleen, Intestine

Available upon request. Contact CIRS customer service at admin@cirsinc.com for more information.

| SLAB SIZE (CM) | THICKNESS (CM) | | | | |
|----------------|----------------|-----|-----|-----|-----|
| | 0.5 | 1.0 | 2.0 | 3.0 | 4.0 |
| 10 X 10 | X | X | X | X | X |
| 20 X 20 | X | X | X | | |
| 30 X 30 | | X | X | X | X |

