

MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007

Fluoroscopy Mesh Phantom

A routine test object designed to provide an ongoing check of the resolution of fluoroscopy systems. Image quality is measured simply by counting the number of mesh patterns resolved in the image. An ongoing record of these numbers will reveal any trend towards deterioration in imaging performance.

Available in three resolution ranges:

• FLUOROMESH 16-60 contains 8x mesh wedges with a range of 16 to 60 lines/inch

FLUOROMESH 30-100

20

40

60

contains 8x mesh wedges with a range of 30 to 100 lines/inch

FLUOROMESH 60-150

contains 8x mesh wedges with a range of 60 to 150 lines/inch





MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007

Fluoroscopy

Product X-rays



fig. 1 FluoroMesh 16-60 x-ray



fig. 2 FluoroMesh 30-100 x-ray



fig. 3 FluoroMesh 60-150 x-ray

深圳市展业达鸿科技有限公司 日先生:15920060912(微信同号) 0755-22934005(座机) 地址:深圳市福田区八卦二路八卦岭工业区615栋419 邮箱:hongqi@thingstet.com Copyright © 2014 Leeds Test Objects Ltd | SpecSh 419 図址; 愛明4 thingstest.com